



## HEINE AllSpec® Disposable Tips

DATA	
Description	HEINE AllSpec® Disposable Tips – single use plastic tips for HEINE Otoscopes, different sizes for infants and adults
Catalogue Number	Infants 2.5 mm dia: Pack of 1000 pcs. B-000.11.128 / Dispenser pack of 250 tips B-000.11.153; Adults 4 mm dia: Pack of 1000 pcs. B-000.11.127 / Dispenser pack of 250 tips B-000.11.152
Version / Date	V01/10.07.2018
GENERAL	
Features	Most secure fit: especially designed for notches of HEINE Otoscopes; best light transmission: matt inside and no stray light due to opaque material; assured quality: no sharp edges; stackable and suitable for Tip dispenser
Dimensions product	Infants: Ø proximal 17.5 mm, Ø distal 2.5 mm, length 36 mm; Adults: Ø proximal 15 mm, Ø distal 4 mm, length 30.5 mm
Weight packaging	Pack of 1000 pcs. 510 g / Dispenser pack of 250 tips 160 g
Material	Polystyrol
REACH/RoHS	Compliant
Phthalate	Product is phthalate free
Latex	Product is latex free
Biocompatibility	Compliant
Surface	Anthracite similar to RAL 7016
Environmental conditions operation	10 °C to 35 °C / 30 % to 75 % rel. hum. / 700 hPa to 1060 hPa
Environmental conditions storage	5 °C to 45 °C / 45 % to 80 % rel. hum. / 500 hPa to 1060 hPa
Environmental conditions transport	-20 °C to 50 °C / 45 % to 80 % rel. hum. / 500 hPa to 1060 hPa
Durability	Single-use in order to prevent cross-contamination
Customs Code	90189084
EAN/GTIN	Pack of 1000 pcs. 2.5 mm dia. (infants) B-000.11.128: 4053755113897; Pack of 1000 pcs. 4 mm dia. (adults) B-000.11.127: 4053755113866
MECHANICAL	
Connections	Fits all HEINE F.O. Otoscopes, with the HEINE Tip Adaptor also suitable for HEINE BETA 100 and HEINE K 100 Otoscopes
Imprints	HEINE, CE-sign molded-in (proximal on the outside)
USE	
General	HEINE Otoscopes are intended for the examination of the ear, the external auditory canal and the tympanum. The HEINE Otoscope may only be inserted into the auditory canal when a tip is fitted.
Instructions for use	Take for each patient and each ear a new tip. Place the HEINE Tip on the otoscope so that the projection inside the tip fits into the bayonet slit. Tighten the tip by twisting it slightly to the right.
Intended user	The instruments may only be used by qualified medical personnel.
CLEANING, DISINFECTION	
Procedure	None, product is single-use

